Bürette / Buret

continuous

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[Operating manual](#)

[Mode d'emploi](#)

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[Istruzione](#)

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**Before using the instrument for the first time,
ensure it is rinsed carefully or discard the first
few samples dispensed.**

VITLAB GMBH

Linus-Pauling-Str. 1
63762 Grossostheim
Germany
T +49 6026 977 99 0
info@vitlab.com
www.vitlab.com

Impressum

Vitlab GmbH

Linus-Pauling-Str. 1
63762 Grossostheim, Germany
tel: +49 6026 97799-0
fax: +49 6026 97799-30
info@vitlab.com
www.vitlab.com

Do you need other operating manuals? Please refer to
<https://www.vitlab.com/manuals> or use the following QR Code:

The original operating manual is written in German. Other languages are translations of the original operating manual.
Technical changes, errors or misprints reserved.

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DECLARATION OF CONFORMITY
– China RoHS 2

VITLAB GMBH has made reasonable efforts to ensure that restricted hazardous substances are not used in VITLAB products above the limits specified by China RoHS 2.

In order to determine the concentration of hazardous substances in all homogeneous materials of the subassemblies, a “Product Conformity Assessment” (PCA) procedure was performed. As defined in GB 26572-2025 the “Maximum Concentration Value” limits (MCV) apply to these restricted substances:

+ Lead (Pb):	0.1 %	+ Hexavalent chromium (Cr(+VI)):	0.1 %
+ Mercury (Hg):	0.1 %	+ Polybrominated biphenyls (PBB):	0.1 %
+ Cadmium (Cd):	0.01 %	+ Polybrominated diphenyl ether (PBDE):	0.1 %
+ Bis(2-ethylhexyl) phthalate (DEHP):	0.1 %	+ Dibutyl phthalate (DBP):	0.1 %
+ Butyl benzyl phthalate (BBP):	0.1 %	+ Diisobutyl phthalate (DIBP):	0.1 %

Environmental Friendly Use Period (EFUP)

EFUP defines the period in years during which the hazardous substances contained in electrical and electronic products will not leak or mutate under normal operating conditions. During normal use by the user such electrical and electronic products will not result in serious environmental pollution, cause serious bodily injury or damage to the user’s assets.

The Environmental Friendly Use Period for VITLAB instruments is 40 years.

Material Content Declaration for VITLAB Products

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Apart from the disclosures in the above table, the subassemblies are not intentionally manufactured or formulated with lead (Pb), mercury (Hg), cadmium (Cd), hexavalent chromium (Cr+VI), polybrominated biphenyls (PBB), polybrominated diphenyl ethers (PBDE), bis(2-ethylhexyl) phthalate (DEHP), butyl benzyl phthalate (BBP), dibutyl phthalate (DBP) and diisobutyl phthalate (DIBP).

Products manufactured by VITLAB may enter into further devices or can be used together with other appliances. With these third party products and appliances in particular, please note the EFUP labeled on these products. VITLAB will not take responsibility for the EFUP of those products and appliances.

Place, date: Grossostheim, 01/07/2026

Peter Schütte
 (Managing Director)

Dr. Stephan Schmidt
 (Regulatory Compliance Manager)

VITLAB GmbH
 Linus-Pauling-Str.1
 63762 Grossostheim
 Germany
 Telefon: +49 6026 97799-0
 Fax: +49 6026 97799-30
 E-mail: info@vitlab.com
 Internet: www.vitlab.com

Contents

1 Introduction	23
1.1 Scope of supply	23
1.2 Using the operating manual	23
2 Safety regulations	23
2.1 General safety regulations	23
2.2 Operating restrictions and function	24
2.3 Application restrictions	25
2.4 Operating exclusions	25
2.5 Recommended Application Range for VITLAB® continuous	26
3 Functional and operational components	27
4 First Steps.....	28
5 Titrating	28
5.1 Preparation for titration	28
5.2 Priming	29
5.3 Titrating.....	29
5.4 Changing the reagent bottle	30
6 Cleaning.....	30
6.1 Emptying	31
6.2 Standard cleaning.....	31
6.3 Intensive Cleaning	31
6.4 Clean Discharge/recirculation valve.....	32
7 Battery replacement	32
8 Checking the Volume (Calibration)	33
9 Recalibration	34
9.1 Recalibration	34
9.2 Reset function	35
10 Technical data	35
11 Troubleshooting.....	36
12 Order Information	36
13 Repair – calibration service.....	37
13.1 Sending for repair.....	37
14 Calibration service.....	37
15 Warranty.....	38
16 Disposal.....	38

1 Introduction

1.1 Scope of supply

Digital Buret VITLAB continuous E / RS with GL 45 connecting thread and threaded adapter made of PP in sizes GL 32, GL 38 and S 40 (buttres thread), telescoping filling tube (200 - 350 mm), telescoping discharge tube (140 - 220 mm), 2 Micro batteries 1,5 V (LR03, AAA), performance certificate and this operating manual.

1.2 Using the operating manual

- Please carefully read the operating manual before using the device for the first time.
- Keep the operating manual in an easily accessible place. It is part of the instrument.
- Be sure to include the operating manual if you transfer possession of this device to a third party.

1.2.1 Signal words and their meaning

⚠ WARNING! ...	The signal word WARNING indicates a dangerous situation that, if not avoided, could result in death or serious injury.
⚠ CAUTION! ...	The signal word CAUTION indicates a hazardous situation that, if not avoided, could result in moderate or minor injury.
NOTICE! ...	The signal word NOTICE is used to address actions that are not related to physical injury. Example: Possible property damage. The note shows operating tricks or other useful information

1.2.2 Presentation of descriptions of actions

Format	Meaning
1. Task	Indicates a task.
a., b., c.	Indicates the individual steps of a task.
>	Indicates a prerequisite for a task.
⇒	Indicates a result of a completed task.

2 Safety regulations

2.1 General safety regulations

Please read carefully!

The can be used in combination with hazardous materials, work processes and equipment. However, the operating manual cannot cover all of the safety issues that may occur in doing so. It is the user's respon-

sibility to ensure compliance with the safety and health regulations and to specify the corresponding restrictions before use.

1. Use the utmost caution when dispensing caustic, poisonous, radioactive or hazardous chemicals.
2. Observe general safety regulations, e. g., wear protective clothing, goggles and gloves.
3. Observe the Operating Manual and information from reagent manufacturers.
4. Never use the instrument in an atmosphere that might be explosive.
5. Use the instrument only for dispensing liquids, with strict regard to the defined Operating Exclusions and Limitations. If in doubt regarding the suitability of the instrument, for a particular application contact the manufacturer.
6. Before use, always verify that the instrument is in good working order, e. g., piston moves smoothly; filling and discharge tube are firmly seated and are properly attached, etc.
7. Never use force on the instrument. Use of force may result in danger to the user or other persons.
8. While dispensing, the discharge tube must always point away from the user or other persons. Avoid splashes. Only dispense into suitable vessel.
9. Never carry the mounted instrument by its upper casing.
10. Clean the instrument before disassembling.
11. Use only original manufacturer's accessories and spare parts. Do not make any technical modifications.
12. In case of difficulty (e. g., hand wheels difficult to move, leakage), immediately stop dispensing. Before further use, repair the instrument as described in this Manual. Contact the manufacturer if necessary.
13. The included 1.5 V micro-batteries are not re-chargeable!

CE labeling

With this mark, we confirm that the product complies with the requirements set out in the EC Directives and has been subjected to the specified testing procedures.

2.2 Operating restrictions and function

The instrument has a continuous, pulse-free dispensing technique and is designed for titrating liquids, observing the following physical limits:

- Working- and charging temperature: +15 °C to + 40 °C (for instrument and liquid)
- Storage temperature: -20 — + 50 °C
- Rel. Humidity 0 % - 90 % non-condensing
- Height above sea level, operation/storage: 2000 m / 12200 m
- Density up to 2.2 g/cm³
- Vapor pressure up to 500 mbar

The dispensed volume of 0,01 ml to 999,9 ml appears on the display.

⚠ WARNING! To avoid splashing do not remove the discharge tube from its support!

2.3 Application restrictions

Liquids, which form deposits may make the piston difficult to move or may cause jamming (e.g., crystallizing solutions or concentrated alkaline solutions).

When dispensing inflammable media, make sure to avoid to buildup of static charge, e.g., do not dispense into plastic vessels; do not wipe instruments with a dry cloth.

⚠ WARNING! If there is a sign of a potential malfunction (e.g., piston difficult to move) never use force. Stop dispensing immediately and perform cleaning, see Cleaning, p. 30. Contact the manufacturer, if necessary.

NOTICE! Compatibility of the instrument for the intended application (i.e., trace analysis) must be carefully checked by the user. Contact the manufacturer, if necessary.

⚠ WARNING! Do not remove the discharge tube from its support! It must not be used as flexible discharge tube!

2.4 Operating exclusions

When the instrument is correctly used, the dispensed liquid comes into contact with only the following materials: Borosilicate glass 3.3, FEP, ETFE, PFA, PTFE and platinum-iridium.

Never use this instrument for:

- Liquids which attack FEP, ETFE, PFA and PTFE
- Solutions containing hydrofluoric acid
- Solutions which tend to crystallize, fuming acids and concentrated bases
- Suspensions containing solid particles
- Solutions which decompose and form solid particles (e.g., Biuret reagent);
- Substances which undergo catalytic transformation or react with Platinum-Iridium (e.g., H_2O_2)
- Carbon disulfide as this media inflames easily
- The instrument should not be used in an aggressive atmosphere (e.g., HCl fumes).
- Autoclaving the instrument is not permitted!



2.5 Recommended Application Range for VITLAB® continuous

The bottletop-buret VITLAB® continuous E/RS can be used for the following titration media (max. conc. 1 mol/l).

Medium	Medium
Ammonium iron (II) sulfate solution	Sodium arsenite solution
Ammonium thiocyanate solution	Sodium carbonate solution
Barium chloride solution	Sodium chloride solution
Bromide bromate solution	Sodium nitrite solution
Cerium (IV) sulfate solution	Sodium thiosulfate solution
EDTA solution	Sodium hydroxide solution*
Iron (II) sulfate solution	Oxalic acid
Acetic acid	Perchloric acid
Potassium hydroxide solution*	Nitric acid*
Potassium bromate solution	Hydrochloric acid*
Potassium bromide bromate solution	Sulphuric acid*
Potassium dichromate solution	Silver nitrate solution*
Potassium iodate solution	Tetra-n-butylammonium hydroxide solution
Potassium permanganate solution*	Zinc sulfate solution
Potassium thiocyanate solution	

* To prevent humidity and carbon dioxide from entering the eluents, the use of an accordingly filled drying tube is recommended. For acids, a molecular sieve is suitable for absorbing air moisture; for alkalis, soda lime is suitable as a weak CO₂ adsorbent.

3 Functional and operational components

4 First Steps

Before the first use please insert the provided batteries:

- a. With your thumb, press on the latch located below the control keys. With your other hand, simultaneously push the front part of the outer casing toward the top.

NOTICE! Do not apply force.

- b. Insert two new micro batteries 1,5 V (LR03/AAA). Ensure that the polarity of the batteries corresponds to the symbols on the battery compartment.
- c. To reassemble the casing, position the bottom edge in place, then gently push downward until you can hear it locking.

5 Titrating

5.1 Preparation for titration

- a. Observe Operating Limitations and general safety regulations.
- b. Adjust the length of the telescopic intake tube (according to the bottle to be used) and push the intake tube into the intake valve as far as it will go. Cut the lower end at an angle.

- c. Screw the instrument (thread GL 45) onto the reagent bottle and align it with the bottle label. For bottles with different thread diameter select a suitable adapter!

⚠ WARNING! Avoid splashing the reagent! Never carry the device by the top casing. Always support both the instrument and the bottle. Carry the mounted instrument only as shown in figure.

NOTICE! If the bottle is shorter than the filling tube, separate the two sections and cut both filling tubes to the correct length.

5.2 Priming

⚠ WARNING! Never point the ejection cannula towards the user or other people. Do not move the hand wheels until the instrument has been correctly mounted, and the valve switch is in the "recirculate" position.

- Set the valve switch to "recirculation."
- Turn the hand wheels forward 5 - 10 turns to remove air from the mechanism.
- Set the valve switch to the normal operating position, and hold a suitable vessel below the discharge tube outlet.
- Turn the hand wheels forward until air bubbles have disappeared in the discharge tube.

5.3 Titrating

- To switch on the LCD, press the start key.
- ⇒ The titration mode appears in the display.
- Place a suitable collecting vessel below the discharge tube outlet.

⚠ WARNING! Wipe the remaining reagent drops from the discharge tube into the vessel.

- Dispense reagent by turning the hand wheels forward until the desired value has been reached. If the hand wheels are rotated backward accidentally, they do not engage and the dispensing process is not affected.

⚠ WARNING! To avoid splashing of reagent, always turn the hand wheels using a slow and steady motion.

- d. To start a new titration, press the "Start" key to reset the display to zero.
- e. The display switches off automatically after approx. 5 minutes of idle time. However, the titrated value remains stored. It will reappear when the "Start" or "Pause" key are pressed, and the titration can be continued.
- f. When the titrating action is finished, set the valve switch back to "recirculation."

NOTICE! When the titrating is completed or interrupted, always set the valve switch to "recirculation" to prevent the accidental release of liquid from the discharge tube.

5.4 Changing the reagent bottle

Changing the reagent bottle during a titration process:

- a. Press the "Pause" key. The titrated value remains stored in the display.
- b. For instructions on changing bottles, see Emptying, p. 31.
- c. Mount a new bottle and prime the instrument, see Priming, p. 29.
- d. Press the "Pause" key again. The "Titration" icon reappears in the display, and titration can be resumed, starting with the stored value.

6 Cleaning

Instrument will only function safely if cleaned regularly. Be sure to clean instrument:

- Immediately, if the hand wheels move stiffly
- Before changing reagents
- Before longer breaks in use
- Before any maintenance and repair work
- Daily after use with solutions prone to crystallization and concentrated bases

⚠ WARNING! The intake mechanism, valves, and intake/discharge tubes are filled with reagent. Observe general safety regulations, e. g., wear protective clothing, goggles and gloves.

6.1 Emptying

- a. Place instrument into a suitable basin.
- b. Unscrew the instrument and lift it high enough so that the intake tube is no longer immersed in liquid.
- c. Cautiously tap the intake tube against the inside of the bottle so that any remaining reagent runs out.
- d. Remove the dispenser from the bottle and mount it on another empty bottle.
- e. Point the discharge tube into the opening of the used bottle, and empty the instrument by turning the hand wheels. Then set the valve switch to "recirculation" and turn the hand wheels again to empty the recirculation channel.

6.2 Standard cleaning

- a. Mount the instrument on a bottle filled with a suitable cleaning solution.
- b. Rinse the instrument thoroughly by turning the hand wheels.
- c. Remove the instrument from the bottle and empty it completely by turning the hand wheels, as described above.
- d. Mount the instrument on a bottle filled with distilled water. Rinse thoroughly, then empty, as described above.

6.3 Intensive Cleaning

Intensive cleaning must be preceded by the standard cleaning, p. 31 procedure. This procedure must be followed if the hand wheels are hard to move, or if the instrument is heavily soiled. For this purpose, the dispenser must be partially disassembled.

⚠ WARNING! Use personal protective equipment.

Before disassemble, always perform the standard cleaning procedure. To prevent injury from chemicals, always wear eye protection, protective clothing and protective gloves. Avoid splashing the reagent.

- a. Pull off the telescopic intake tube and clean it with a soft bottle brush. Replace if necessary, especially if cracked or worn.
- b. Clean discharge tube with a soft brush. Replace if necessary, especially if cracked or worn.

-
-
- c. To clean the discharge incl. recirculation valve, see Clean Discharge/recirculation valve, p. 32
 - d. **NOTICE!** The intake mechanism including the internal valves can only be serviced by trained and authorized service personnel. If continues to be soiled after repeated rinsing, please send the instrument to your authorized dealer for repair, see Repair – calibration service, p. 37.

6.4 Clean Discharge/recirculation valve

⚠ WARNING! Never use force during disassembly and assembly. Verify that all components fit tightly and securely.

- a. Set valve switch to „recirculation“.
- b. Remove the valve switch, then lift the discharge tube support.
- c. Unscrew the discharge tube securing nut and pull out the discharge tube.
- d. Pull out the discharge/recirculation valve.
- e. Clean all parts. Replace if necessary.
- f. To reassemble, push in the discharge/recirculation valve by hand.
- g. Push in the discharge tube as far as possible and tighten the securing nut.
- h. Mount the discharge tube support (16) and the valve switch.

NOTICE! Subsequent to reassembly, gravimetrically check volume see Checking the Volume (Calibration), p. 33.

7 Battery replacement

When battery power runs low, a crossed-out battery symbol appears at the top left of the display. The batteries must then be replaced.

- a. With your thumb, press on the latch (7) located below the control keys. With your other hand, simultaneously push the front part of the outer casing toward the top.

NOTICE! Never use force!

- b. Remove the old batteries and dispose of properly.
- c. Insert two new micro batteries 1,5 V (LR03/AAA). Ensure that the polarity of the batteries corresponds to the symbols on the battery compartment.

-
-
- d. To reassemble the casing, position the bottom edge in place, then gently push downward until you can hear it locking.

8 Checking the Volume (Calibration)

Depending on use, we recommend that gravimetric testing of the instrument be carried out every 3-12 months. Set the 10 ml device to 3 decimal places for calibration. This time frame should be adjusted to correspond with individual requirements. The complete testing procedure (SOP) can be downloaded at www.brand.de. You can also carry out a functional test at more frequent intervals, e.g., by titrating against a standard solution.

Gravimetric volume testing according to DIN EN ISO 8655- 6 (for measurement conditions, see Fehlergrenzen) is performed as follows:

1. Prepare the instrument

Clean the device (Cleaning), fill with test fluid (distilled water) and air carefully.

2. Check the volume

- a. Dispense 5 drops in a separate vessel and wipe off the titration tube tip.
- b. Press the CLEAR button to set the display value to zero.
- c. At 10 pipetting series and weighings in 3 volume ranges (100 %, 50 %, 10 %) are recommended.
- d. Turn the hand wheels with both hands without stopping until the test volume is shown on the display. Wipe off the tip of titration tube.
- e. Weigh the pipetted amount with an analysis scale. (Please refer to the user manual of the scale manufacturer.)
- f. Calculate the dispensed volume. The Z factor takes account of the temperature and air buoyancy.

Calculation (for nominal volume)

x_i = Weighing results

N = Number of weighings

V_0 = Nominal volume

Z = Correction factor (e.g. 1.0029 µl/mg at 20 °C, 1013 hPa)

Mean value:

Mean volume:

Accuracy*:

$$\bar{x} = \frac{\sum x_i}{n}$$

$$\bar{V} = \bar{x} * Z$$

$$A\% = \frac{\bar{V} - V_0}{V_0} * 100$$

Variation coefficient*:

$$CV\% = \frac{100 s}{\bar{V}}$$

Standard deviation*:

$$s = Z * \sqrt{\frac{\sum (x_i - \bar{x})^2}{n - 1}}$$

*) = Calculation for accuracy (A%) coefficient of variation (CV%): A% and CV% are calculated using the formulas of statistical quality control.


NOTICE! Test instructions (SOPs) are available for download at www.vitlab.de.

9 Recalibration

9.1 Recalibration

After prolonged use, or for specific applications instrument can be adjusted in order to compensate for differences in accuracy.

- After the titration volume has been dispensed (e. g., 25 ml to 50 ml), press the Pause key. The last volume which has been titrated is then stored.

NOTICE! It is not possible to perform calibration with a value < 10 ml or > 90 ml, as the test volume is then too low or too high. At the attempt to enter the CAL mode with a volume < 10 ml or > 90 ml, the following displays will appear:

- Hold down the ▲-key and the ▼-key simultaneously for 3 seconds until "CAL" starts flashing in the top right-hand corner of the display.
- Using the ▲ or ▼ key, change the value shown in the display to the value of the mean measured volume, see Checking the Volume (Calibration), p. 33.
- Press the "Start" key.

⇒ The adjusted value is automatically accepted; the display is reset to zero, and the calibration procedure is completed. The symbol "C" appears in the display from now

on, indicating that the factory calibration has been modified. (To revert to the factory calibration, see Reset function, p. 35)

NOTICE! When in titrating mode and the "C" symbol is displayed, the current calibration value can be checked by pressing the ▲ or ▼ key. When the keys are released, the display automatically reverts back to the last value displayed.

9.2 Reset function

The reset function restores the factory calibration.

- Press the "Start" key. The display is reset to zero, and the instrument changes into titrating mode.
 - Press the "Start" key and the ▼ key simultaneously for 3 s.
- ⇒ The "C" symbol disappears from the display, and the original factory calibration is restored.

10 Technical data

	Volume per turn**	Systematic error*, A%	Random error*, CV%
continuous E	2.5 ml	$\leq \pm 0.2$	≤ 0.1
continuous RS	5 ml	$\leq \pm 0.2$	≤ 0.1

* Values based on nominal volume

** Dispensing volume per rotation of the hand wheels

Adjustment conditions: Distilled water. EX 20 ± 0.5 °C

Number of testing procedures: 10 according to DIN EN ISO 8655

(Technical changes reserved!)

Power and current data

- Operating voltage 3V (2 LR03/AAA, 1.5 V batteries each)
- Power consumption less than 10 mA

11 Troubleshooting

Fault	Possible causes	Corrective action
Air is sucked in or air bubbles in the discharge cannula.	The instrument has not been adequately vented.	Follow the procedure for bleeding the Instrument, see Priming, p. 29.
	The suction pipe is not sufficiently pushed or damaged.	Carry out the cleaning procedure, see Cleaning, p. 30. Push the suction pipe all the way onto the suction valve, shorten it if necessary or replace it if necessary.
	The end of the aspiration pipe does not immerse itself in the fluid.	Extend the suction pipe until it is immersed in the liquid.
No fluid is aspirated.	The suction valve is clogged or stuck.	For cleaning procedures, see Cleaning, p. 30. If applicable. Return the unit for repair, see Repair – calibration service, p. 37.
The titration volume delivered is too low.	The suction valve is dirty or damaged.	Perform intensive cleaning, see Intensive Cleaning, p. 31. If necessary return the instrument for repair, see Repair – calibration service, p. 37.
	The suction pipe is not sufficiently pushed or damaged.	Push the suction pipe all the way onto the suction valve, shorten it if necessary or replace it if necessary.
	The buret is in "C" mode.	Perform a reset to return to the factory setting, see Reset function, p. 35.

If these measures do not help, return instrument for repair.

12 Order Information

	Cat. No.
continuous E	1620506
continuous RS	1620507

Accessories and replacement parts

Article	Cat. No.
Telescopic intake tube (200 - 350 mm) (FEP, ETFE, PTFE)	1671085
Discharge/recirculation valve (PTFE, PFA) for continuous E und RS	1655085
Dosing cannula, complete	1650135
Discharge tube support, adjustable	1650162
Dry tube, complete (without filling)	1671095
Micro batteries 1.5 V, (LR03/AAA), 2 pieces	1670216
Plastic tripod, PP	1671116

Threaded bottles, coated

Volume	Thread	Part No.
1000 ml	GL 45	1671500
2500 ml	GL 45	1671510

Threaded adapter

Thread	Part No.
GL 45- GL 32	1670180
GL 45- GL 38	1670110
GL 45- S 40	1670120
GL 32 - NS 19/26	1670066
GL 32 - NS 24/29	1670067
GL 32 - NS 29/32	1670068

13 Repair – calibration service

If any instrument malfunction cannot be remedied in the lab by simply exchanging replacement parts, the instrument must be sent in for repair.

Please note that for safety reasons only clean and decontaminated instruments can be inspected and repaired!

13.1 Sending for repair

NOTICE! Transporting hazardous materials without approval is prohibited by law.

1. Clean the instrument thoroughly and decontaminate.
2. Complete the “Declarations of Health Safety Compliance” form (forms can be requested from the dealer or manufacturer or downloaded from www.vitlab.com).
3. Send the completed form together with the instrument to the manufacturer or dealer. Include an exact description of the malfunction and the media used.

Shipment is at the risk and the cost of the sender.

14 Calibration service

The ISO 9001 and GLP guidelines require regular inspection of your volume measuring instruments. We recommend performing a volume check every 3 to 12 months. The cycle is dependent on the individual requirements of the instrument. Checks should be performed more frequently, in case of high frequency of use or the use of aggressive media.

The detailed testing procedure is available for download at www.vitlab.com.

VITLAB also offers the option to have your instruments calibrated through our calibration service.

Simply send us the instruments to be calibrated, accompanied by details about which type of calibration you would like. The instruments will be returned to you after a few days together with a test report (factory calibration) or a DAkkS calibration certificate. More information can be obtained from your dealer or directly from VITLAB.

15 Warranty

We shall not be liable for the consequences of improper handling, use, servicing, operating or unauthorized repairs of the device or for the consequences of normal wear and tear, especially of wearing parts such as pistons, seals, valves and the breakage of glass. The same applies for failure to follow the instructions of the operating manual. We are not liable for damage resulting from disassembly beyond that described in the operating manual or if non-original spare parts or components have been installed.

16 Disposal

The icons   mean: Electrical appliances as well as batteries and rechargeable batteries must not be disposed of in household waste (unsorted municipal waste). Dispose of these separately, e.g. at a collection point.

⚠ WARNING! Risk of fire due to heating after short circuit:

Do not short-circuit batteries or rechargeable batteries to discharge! Tape the plug with adhesive tape to prevent short circuits in the collection container! Never disassemble the batteries!

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