Micro dispenser, fix | Micro dispenser, dual volume

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Before using the instrument for the first time, ensure it is rinsed carefully or discard the first few samples dispensed.

VITLAB GMBH

Linus-Pauling-Str. 1
63762 Grossostheim
Germany
T +49 6026 977 99 0
info@vitlab.com
www.vitlab.com

Impressum

Vitlab GmbH

Linus-Pauling-Str. 1
63762 Grossostheim, Germany
tel: +49 6026 97799-0
fax: +49 6026 97799-30
info@vitlab.com
www.vitlab.com

Do you need other operating manuals? Please refer to
<https://www.vitlab.com/manuals> or use the following QR Code:

The original operating manual is written in German. Other languages are translations of the original operating manual.
Technical changes, errors or misprints reserved.

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Contents

1	Introduction	19
1.1	Scope of supply	19
1.2	Using the operating manual	19
2	Safety regulations	19
2.1	General safety regulations	19
2.2	Target group	20
2.3	Improper use	20
2.4	Foreseeable misuse	21
2.5	Function and usage limits	21
2.6	Operating exclusions	21
2.7	Application restrictions	22
2.8	Storage conditions	22
3	Functional and operational components	23
4	Dispensing	24
4.1	Preparation for dispensing	24
4.2	Priming	24
4.3	Switch volume (only piccolo 2) .	24
4.4	Dispensing	24
5	Cleaning	25
5.1	Emptying	25
5.2	Standard cleaning	25
5.3	Intensive cleaning	26
5.4	Cleaning/replacing the valves ...	26
6	Checking the Volume (Calibration)	27
7	Troubleshooting	28
8	Product markings	29
9	Technical data	29
10	Order Information	30
11	Repair – calibration service	30
11.1	Sending for repair	30
12	Calibration service	31
13	Warranty	31
14	Disposal	31

1 Introduction

1.1 Scope of supply

VITLAB® piccolo 1 or VITLAB® piccolo 2 complete with discharge tube, intake tube and tool.

1.2 Using the operating manual

- Please carefully read the operating manual before using the device for the first time.
- Keep the operating manual in an easily accessible place. It is part of the instrument.
- Be sure to include the operating manual if you transfer possession of this device to a third party.

1.2.1 Signal words and their meaning

DANGER! ...	The signal word DANGER indicates a dangerous situation that, if not avoided, will result in death or serious injury.
⚠ WARNING! ...	The signal word WARNING indicates a dangerous situation that, if not avoided, could result in death or serious injury.
⚠ CAUTION! ...	The signal word CAUTION indicates a hazardous situation that, if not avoided, could result in moderate or minor injury.
NOTICE! ...	The signal word NOTICE is used to address actions that are not related to physical injury. Example: Possible property damage. The note shows operating tricks or other useful information

1.2.2 Presentation of descriptions of actions

Format	Meaning
1. Task	Indicates a task.
a., b., c.	Indicates the individual steps of a task.
>	Indicates a prerequisite for a task.
⇒	Indicates a result of a completed task.

2 Safety regulations

2.1 General safety regulations

Please read carefully!

The can be used in combination with hazardous materials, work processes and equipment. However, the operating manual cannot cover all of the safety issues that may occur in doing so. It is the user's respon-

sibility to ensure compliance with the safety and health regulations and to specify the corresponding restrictions before use.

- Every user must read and observe this operating manual before using the instrument.
- Follow the general hazard instructions and safety regulations, e.g. wear protective clothing, eye protection and protective gloves.
- Follow the instructions given by the reagent manufacturer.
- Use the instrument only for dispensing liquids, with strict regard to the defined limitations of use and operating limitations.
- Comply with the operating exclusions; see Operating exclusions, p. 21. In case of doubt, contact the manufacturer or dealer.
- Never use force.
- Always perform work in a manner that does not endanger users or other people. When dispensing, the discharge tube must always point away from you or any other person. Avoid splattering. Use only suitable vessels.
- Clean the instrument before disassembling.
- Use only original manufacturer's accessories and spare parts. Do not make any technical modifications. Do not dismantle the instrument any further than is described in the operating manual!
- Always check that the device is in proper working condition before use. If there is a sign of a potential malfunction (i.e., piston difficult to move, sticking valves or leakage), immediately stop dispensing and consult Troubleshooting, p. 28. Contact the manufacturer, if necessary.

2.2 Target group

The operating manual is intended for users who use the laboratory instrument in the course of their professional activities. Users are familiar with the typical safety regulations and working methods in laboratories and have been trained accordingly. They can recognize potential hazards and protect themselves from them. The operating manual assumes this expertise and does not replace basic laboratory training or specific safety training.

2.3 Improper use

Various risks may arise if the laboratory instrument is used improperly. These risks include: inaccurate liquid delivery, damage to the laboratory instrument, and the risk of contamination, infection, and injury from contact with the dispensed media.

Any use of the laboratory instrument other than for dispensing liquids within the defined operating limits is considered improper use.

2.4 Foreseeable misuse

A typical misuse is insufficient cleaning with crystallizing or hygroscopic liquids.

2.5 Function and usage limits

This instrument is designed for dispensing of aqueous and highly diluted agents, observing the following physical limits:

- Useage temperature from +15 °C to +40 °C of instrument and reagent
- Density up to 2.2 g/cm³
- Vapor pressure up to 500 mbar
- Viscosity up to 200 mm²/s

2.6 Operating exclusions

NOTICE! The user is responsible for checking the compatibility of the instrument with the intended application.

When the instrument is correctly used, the dispensed liquid comes into contact with only the following materials:

- Borosilicate glass, FEP, ETFE, PFA, PTFE and Platinum-Iridium.

Never use the dispenser for:

- Concentrated acids and alkalis
- Liquids which attack FEP, ETFE, PFA and PTF; e.g., dissolved sodium azide
- Liquids that attack borosilicate glass (e.g., hydrofluoric acid)
- Suspensions (e.g., activated charcoal) because solid particles may clog or damage the instrument
- Solutions which decompose and form solid particles (e.g., Biuret reagent);
- Substances which undergo catalytic transformation or react with Platinum-Iridium (e.g., H₂O₂)
- Explosive liquids (e.g., carbon disulfide)
- Tetrahydrofuran
- Trifluoroacetic acid

2.7 Application restrictions

Liquids, which form deposits may make the piston difficult to move or may cause jamming (e.g., crystallizing solutions or concentrated alkaline solutions). If the piston movement becomes sluggish or stiff, the instrument should be cleaned immediately (Cleaning, p. 25).

When dispensing inflammable media, make sure to avoid to buildup of static charge, e.g., do not dispense into plastic vessels; do not wipe instruments with a dry cloth.

The instrument is designed for general laboratory applications and complies with the requirements of the relevant standards (e.g., DIN EN ISO 8655). Compatibility of the instrument for a specific application (e.g., trace material analysis or food sector) must be checked by the user. Approvals for specific applications, e.g. for production and administration of food, pharmaceuticals or cosmetics are not available.


2.8 Storage conditions

Store the unit and accessories in a cool and dry place in cleaned condition only.

Storage temperature from -20 to 50°C (-4 to 122°F).

3 Functional and operational components

Figure: piccolo 1

4 Dispensing


4.1 Preparation for dispensing

- Observe Operating Limitations and general safety regulations.
- Push the intake tube into the intake valve as far as possible and cut the lower end at an angle.
- Screw the dispenser on to the reagent bottle and align it with the bottle label.
⚠ WARNING! Avoid splashing the reagent!
- Place a suitable vessel below the discharge tube.
- The discharge tube must point away from the user at all times.

4.2 Priming

- Place a suitable vessel below the discharge tube.
⚠ WARNING! Reagent may splash!
- Gently push down the piston several times until air bubbles have disappeared from the discharge tube.

4.3 Switch volume (only piccolo 2)

- To set the volume, push down the operating knob and twist handle 90°.
- ⇒ The nose of the housing indicates the volume after releasing operating knob.

4.4 Dispensing

- Place a suitable vessel below the discharge tube.

-
-
- b. Slowly and evenly push the piston down to the stop, then slowly and evenly let the piston slide back up to the upper stop. No air bubbles should be visible during dispensing.

⚠ WARNING! If the piston moves stiffly or is difficult to move, stop dispensing immediately and clean the instrument, see Cleaning, p. 25.

⚠ WARNING! Reagent may escape from the discharge tube.

5 Cleaning

Instrument will only function safely if cleaned regularly. Be sure to clean the instrument:

1. Immediately, if the piston moves stiffly;
2. Before changing reagents
3. Before longer breaks in use
4. Before any maintenance and repair work

⚠ WARNING! Danger from media The device, intake tube and ejection cannula may be filled with reagents. To avoid injuries from chemicals:

Wear eye protection, protective clothing and protective gloves.

Always point the discharge tube outlet away from the user or other people.

5.1 Emptying

- a. While the dispenser is still mounted on the bottle, place into a suitable basin.
- b. Unscrew the dispenser and lift it high enough so that the intake tube is no longer immersed in liquid.
- c. Cautiously tap the intake tube against the inside of the bottle so that any remaining reagent runs out.
- d. Remove the dispenser from the bottle and mount it on another empty bottle.
- e. Hold the discharge tube outlet over the opening of the first bottle. Empty the discharge tube by repeated dispensing movements.

5.2 Standard cleaning

- a. Mount the empty dispenser on a bottle filled with a suitable cleaning solution.

⚠ WARNING! Reagent may splash!

- b. Rinse the dispenser by repeated dispensing movements. Place a suitable vessel below the discharge tube.

-
-
- c. Empty the device completely by dispensing several doses, see Emptying, p. 25.
 - d. Screw the device onto a bottle filled with distilled water, rinse thoroughly and then empty, see Emptying, p. 25.

5.3 Intensive cleaning

Intensive cleaning must be preceded by a standard cleaning. This procedure must be followed if the piston is difficult to move. For this purpose, the dispenser must be partially disassembled.

⚠ WARNING! Before disassembling, always perform the standard cleaning procedure. To prevent injury from chemicals, always wear eye protection, protective clothing and protective gloves. Avoid splashing the reagent.

- a. Pull off the intake tube and clean it with a soft bottle brush. Replace if necessary, especially if cracked or worn.
- b. Use a small screw driver to gently pry between the securing ring slot and valve head and carefully push up the securing ring.
- c. Disassemble outer sleeve by pressing down operating knob and pulling the securing ring simultaneously.
- d. Remove the operating knob with the piston and the lifting spring.
- e. Clean the piston and the cylinder with a soft bottle brush and rinse with water.

NOTICE! Insert the piston at a slight angle to avoid damage to the lip seal.

- f. Mount the lifting spring and carefully reinsert the piston.
- g. Replace the outer sleeve and press it down until it audibly locks in place.
- h. Push down the securing ring.

5.4 Cleaning/replacing the valves

5.4.1 Cleaning/replacing the intake valve

NOTICE! Valve ball may fall out when valve is removed

- a. Unscrew the intake valve using the valve tool.
- b. Clean or replace the valve.

-
-
- c. Make sure that valve ball is still in place. To reassemble, handtighten valve and complete using valve tool.

5.4.2 Cleaning/replacing the discharge valve

- a. Remove the discharge tube support. Pull out securing clip with a small screwdriver.
- b. Pull out the discharge valve.
- c. Unscrew the discharge tube securing nut and pull out the discharge tube.
- d. Clean or replace the valve.
- e. To reassemble push in valve by hand and insert securing clip.
- f. Push in the discharge tube as far as possible and tighten the securing nut.
- g. In reverse order mount the discharge tube support.

⚠ WARNING! Never use force during disassembly and assembly. Verify that all components fit tightly and securely. Subsequent to reassembly, gravimetrically check volume see Checking the Volume (Calibration), p. 27.

6 Checking the Volume (Calibration)

Depending on use, we recommend that gravimetric testing of the instrument be carried out every 3–12 months. This time frame should be adjusted to correspond with individual requirements. Gravimetric volume testing according to DIN EN ISO 8655- 6 (for measurement conditions, see Error limits) is performed as follows:

1. Prepare the instrument

Clean the device (Cleaning), fill with test fluid (distilled water) and air carefully.

2. Check the volume

- a. Execute 10 doses with test fluid in 3 volume ranges (100%, 50%, 10%).
- b. For discharge, depress the piston slowly and steadily without force until the lower stop
- c. Wipe off the tip of titration tube.
- d. Weigh the pipetted amount with an analysis scale. Please refer to the manufacturer's instructions for use.
- e. Calculate the dispensed volume. The Z factor takes account of the temperature and air buoyancy.

Calculation (for nominal volume)

x_i = weighing results

n = number of weighings

V_0 = nominal volume

Z = Correction factor (e.g. 1.0029 $\mu\text{L}/\text{mg}$ at 20°C, 1013 hPa)

Mean:

$$\bar{x} = \frac{\sum x_i}{n}$$

Mean volume:

$$\bar{V} = \bar{x} * Z$$

Accuracy*:

$$A\% = \frac{\bar{V} - V_0}{V_0} * 100$$

Standard deviation*:

$$s = Z * \sqrt{\frac{\sum (x_i - \bar{x})^2}{n - 1}}$$

Coefficient of variation*:

$$CV\% = \frac{100 s}{\bar{V}}$$

*) Accuracy and coefficient of variation are calculated according to the formulas of statistical quality control.

NOTICE! Test instructions (SOPs) are available for download at www.vitlab.de.

7 Troubleshooting

Fault	Possible causes	Corrective action
Piston sluggish	Crystal deposits	Clean (see Cleaning, p. 25)
Reagent is not aspirated	Filling valve stuck	Cleaning/replacing the valves, see Cleaning/replacing the valves, p. 26
Dispensing not possible	Discharge valve damaged	Cleaning/replacing Discharge valve valves, see Cleaning/replacing the valves, p. 26
Intake-/discharge tube not firmly connected	Tubes damaged	Intake and Discharge tube, see Cleaning/replacing the valves, p. 26
Air is aspirated	Intake valve not firmly connected	Cleaning/replacing the valves, see Cleaning/replacing the valves, p. 26
	Intake tube not firmly connected or damaged	Push intake tube on firm-ly; replace if necessary
	Reagent with high vapor pressure has been drawn in too quickly	Aspirate reagent more slowly
Dispensed volume too low	Intake valve leaks	Cleaning/replacing the valves, see Cleaning/replacing the valves, p. 26
	Intake tube not firmly connected or damaged	Push intake tube on firm-ly; replace if necessary

If these measures do not help, return instrument for repair.

8 Product markings

Symbol or number	Meaning
	General warning symbol
	Follow the operating instructions
	Use eye protection
	Use hand protection
	Wear protective clothing
www.vitlab.com/ip	Patent information

9 Technical data

VITLAB® piccolo 1 (micro dispenser, fix), GL 28

Volume	Systematic error* A%	Random measurement deviation* CV%
100 µl	≤ ± 3.0	≤ 0.4
200 µl	≤ ± 2.0	≤ 0.4
250 µl	≤ ± 2.0	≤ 0.4
500 µl	≤ ± 1.5	≤ 0.3
1000 µl	≤ ± 1.0	≤ 0.2

* Values based on nominal volume.

VITLAB® piccolo 2 (micro dispenser, variable), GL 28

Volume	Systematic error* A%	Random measurement deviation* CV%
100/ 250 µl	≤ ± 2.0	≤ 0.4
500/ 1000 µl	≤ ± 1.0	≤ 0.2
1000/ 2000 µl	≤ ± 1.0	≤ 0.2

* Values based on nominal volume.

Adjustment conditions: Distilled water. EX 20 ± 0.5 °C

Number of testing procedures: 10 according to DIN EN ISO 8655

(Technical changes reserved!)

10 Order Information

VITLAB® piccolo 1 (micro dispenser, fix), GL 28

Volume	Cat. No.
100 µl	1610501
200 µl	1610502
250 µl	1610503
500 µl	1610504
1000 µl	1610506

VITLAB® piccolo 2 (micro dispenser, variable), GL 28

Volume	Cat. No.
100/ 250 µl	1611503
500/ 1000 µl	1611506
1000/ 2000 µl	1611509

Accessories

Accessories	Cat. No.
Intake tube 150 mm long	1650010
Intake tube 180 mm long	1650011
Intake valve with ball (ETFE/borosilicate glass)	1655011
Discharge tube	1650145
Ejection valve (PFA/borosilicate glass/ platinum-iridium)	1655090
Threaded bottle, 100 ml, coated, GL28	1671505

11 Repair – calibration service

If any instrument malfunction cannot be remedied in the lab by simply exchanging replacement parts, the instrument must be sent in for repair.

Please note that for safety reasons only clean and decontaminated instruments can be inspected and re-paired!

11.1 Sending for repair

NOTICE! Transporting hazardous materials without approval is prohibited by law.

1. Clean the instrument thoroughly and decontaminate.
2. Complete the “Declarations of Health Safety Compliance” form (forms can be requested from the dealer or manufacturer or downloaded from www.vitlab.com).

-
-
3. Send the completed form together with the instrument to the manufacturer or dealer. Include an exact description of the malfunction and the media used.

Shipment is at the risk and the cost of the sender.

12 Calibration service

The ISO 9001 and GLP guidelines require regular inspection of your volume measuring instruments. We recommend performing a volume check every 3 to 12 months. The cycle is dependent on the individual requirements of the instrument. Checks should be performed more frequently, in case of high frequency of use or the use of aggressive media.

The detailed testing procedure is available for download at www.vitlab.com.

VITLAB also offers the option to have your instruments calibrated through our calibration service.

Simply send us the instruments to be calibrated, accompanied by details about which type of calibration you would like. The instruments will be returned to you after a few days together with a test report (factory calibration) or a DAkkS calibration certificate. More information can be obtained from your dealer or directly from VITLAB.

13 Warranty

We shall not be liable for the consequences of improper handling, use, servicing, operating or unauthorized repairs of the device or for the consequences of normal wear and tear, especially of wearing parts such as pistons, seals, valves and the breakage of glass. The same applies for failure to follow the instructions of the operating manual. We are not liable for damage resulting from disassembly beyond that described in the operating manual or if non-original spare parts or components have been installed.

14 Disposal

At the end of their service life, ensure the instruments are properly disposed of and, where possible, recycled. Disposal must be carried out in accordance with applicable national and local regulations.

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