Pipettierhelfer | Pipette controllers

pipeo[®]

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[Operating manual](#)

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[Mode d'emploi](#)

[Istruzione](#)

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DECLARATION OF CONFORMITY
– China RoHS 2

VITLAB GMBH has made reasonable efforts to ensure that restricted hazardous substances are not used in VITLAB products above the limits specified by China RoHS 2.

In order to determine the concentration of hazardous substances in all homogeneous materials of the subassemblies, a “Product Conformity Assessment” (PCA) procedure was performed. As defined in GB 26572-2025 the “Maximum Concentration Value” limits (MCV) apply to these restricted substances:

+ Lead (Pb):	0.1 %	+ Hexavalent chromium (Cr(+VI)):	0.1 %
+ Mercury (Hg):	0.1 %	+ Polybrominated biphenyls (PBB):	0.1 %
+ Cadmium (Cd):	0.01 %	+ Polybrominated diphenyl ether (PBDE):	0.1 %
+ Bis(2-ethylhexyl) phthalate (DEHP):	0.1 %	+ Dibutyl phthalate (DBP):	0.1 %
+ Butyl benzyl phthalate (BBP):	0.1 %	+ Diisobutyl phthalate (DIBP):	0.1 %

Environmental Friendly Use Period (EFUP)

EFUP defines the period in years during which the hazardous substances contained in electrical and electronic products will not leak or mutate under normal operating conditions. During normal use by the user such electrical and electronic products will not result in serious environmental pollution, cause serious bodily injury or damage to the user’s assets.

The Environmental Friendly Use Period for VITLAB instruments is 40 years.

Material Content Declaration for VITLAB Products

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Apart from the disclosures in the above table, the subassemblies are not intentionally manufactured or formulated with lead (Pb), mercury (Hg), cadmium (Cd), hexavalent chromium (Cr+VI), polybrominated biphenyls (PBB), polybrominated diphenyl ethers (PBDE), bis(2-ethylhexyl) phthalate (DEHP), butyl benzyl phthalate (BBP), dibutyl phthalate (DBP) and diisobutyl phthalate (DIBP).

Products manufactured by VITLAB may enter into further devices or can be used together with other appliances. With these third party products and appliances in particular, please note the EFUP labeled on these products. VITLAB will not take responsibility for the EFUP of those products and appliances.

Place, date: Grossostheim, 12/11/2025

W. Nicolaus

Wolfgang Nicolaus
(Managing Director)

D. Schmidt

Dr. Stephan Schmidt
(Compliance Manager)

VITLAB GmbH
Linus-Pauling-Str.1
63762 Grossostheim
Germany
Telefon: +49 6026 97799-0
Fax: +49 6026 97799-30
E-mail: info@vitlab.com
Internet: www.vitlab.com

Table of contents

1 Using the operating manual	20
1.1 Signal words and their meaning..	20
1.2 Presentation of descriptions of actions	20
2 Safety.....	20
2.1 Function.....	22
2.2 Limitations of use	22
2.3 Operating exclusions	22
2.4 Storage conditions.....	22
3 Functional and operating elements.....	23
4 First steps - Connecting the battery.....	23
5 Pipetting.....	24
6 Charging the battery	26
6.1 Replacing the battery	27
7 Replace filter and clean adapter.....	28
8 Cleaning UV disinfection Autoclaving	30
9 Troubleshooting.....	30
10 Ordering Information	31
10.1 Accessories and replacement parts.....	31
11 Warranty.....	31
12 Repair – calibration service.....	32
12.1 Sending for repair.....	32
13 Disposal.....	32

1 Using the operating manual

- Please carefully read the operating manual before using the device for the first time.
- Keep the operating manual in an easily accessible place. It is part of the instrument.
- Be sure to include the operating manual if you transfer possession of this device to a third party.
- You can find updated versions of the Instructions for Use at <https://www.vitlab.com/manuals>.

1.1 Signal words and their meaning

⚠ WARNING or ⚠ WARNING! ...	The signal word WARNING indicates a hazardous situation that, if not avoided, could result in death or serious injury.
⚠ CAUTION or ⚠ CAUTION! ...	The signal word CAUTION indicates a hazardous situation that, if not avoided, could result in minor or moderate injury.
NOTICE or NOTICE! ...	The signal word NOTICE is used to address practices not related to physical injury. Example: Possible property damage.

1.2 Presentation of descriptions of actions

Format	Meaning
1. Task	Indicates a task.
a., b., c.	Indicates the individual steps of a task.
>	Indicates a prerequisite for a task.
⇒	Indicates a result of a completed task.

2 Safety

Please read carefully!

The can be used in combination with hazardous materials, work processes and equipment. However, the operating manual cannot cover all of the safety issues that may occur in doing so. It is the user's responsibility to ensure compliance with the safety and health regulations and to specify the corresponding restrictions before use.

- Every user must read and observe the operating manual enclosed with this laboratory device before and during use. The device may only be used by trained and qualified personnel.
- Follow the general hazard instructions and safety regulations (e.g. wear protective clothing, eye protection and protective gloves).
- When working with infectious or hazardous samples/media (e.g. hazardous substances), the general safety rules in the laboratory must be observed and regulations for handling the samples/me-

dia must be followed. The information provided by the media manufacturer (e.g. safety data sheets) must be observed.

- Use the laboratory device only for pipetting or dispensing media within the defined limitations and restrictions of use. Comply with the operating exclusions.
- If working with inflammable media, make sure to take the necessary precautions to avoid the buildup of static electricity (e.g., do not dispense into plastic vessels and do not wipe instruments with a dry cloth). Do not use the device in potentially explosive atmospheres. If in doubt, contact the manufacturer or supplier.
- Always check that the instrument is in proper working condition before use. If device malfunctions are indicated (e.g. sluggish pistons, leaks or power supply faults), stop working with the device immediately and refer to the “Troubleshooting” section in the operating manual. Contact the manufacturer, if necessary.
- Always perform work in a manner that does not endanger users or other people. Avoid splashes. Only use suitable vessels. Never use unnecessary force when operating, cleaning or servicing the laboratory device.
- If the laboratory device is supplied with power through a power adapter, batteries or accumulators, the components and connections on the device must be regularly checked for proper working condition. Do not operate the laboratory device and its accessories in an unprotected, humid or wet environment.
- Do not make any technical modifications. Use only original replacement parts of the manufacturer; even power adapters or accumulators of identical size and specification from other manufacturers are not permitted. Do not disassemble the laboratory device and its accessories (e.g., power adapters, cables, racks, accumulators or batteries) any further than described in the operating manual!
- Only autoclave the laboratory device when this is permissible according to the operating manual.
- Use only the original power adapter to charge the nickel-metal hydride battery.
- Protect the power adapter from moisture and only use it in conjunction with this device.
- Stop using the power adapter if the connection cable is damaged.
- Only authorized service technicians are permitted to carry out repairs and internal maintenance on the instrument.

 WARNING

Potential risk of explosion due to damaged battery

In extreme cases, improper handling of the instrument or the battery (short circuit, mechanical damage, overheating, etc.) may cause the battery to explode.

2.1 Function

The device serves as an aid for pipetting liquids with graduated and volumetric pipettes made of glass or plastic in the volume range from 0.1 ml to 200 ml and with a suction tube outer diameter < 9.2 mm. When the device is handled correctly, the pipetted liquid comes into contact with the pipette only.

2.2 Limitations of use

This device is intended for pipetting liquids within the following physical limitations:

- from 10 to 40°C (from 50 to 104°F) (device and reagent)
- vapor pressure up to 500 mbar. Absorb slowly above 300 mbar in order to prevent the liquid from boiling.
- Density up to 9 g/cm³

2.3 Operating exclusions

The instrument must not be used for liquids with vapors that can cause corrosion or attack the materials silicone or EPDM.

The instrument is not suitable for use with Pasteur pipettes.

⚠ WARNING	
	Operating the instrument in a potentially explosive atmosphere The instrument must never be operated or charged in a potentially explosive atmosphere. Pipetting highly flammable media (flash point under 0 °C (32 °F), e.g. ether, acetone) is prohibited.


2.4 Storage conditions

Store the unit and accessories in a cool and dry place in cleaned condition only.

Storage temperature from -20 to 50°C (-4 to 122°F).

Relative humidity: 5 % to 95 %.


3 Functional and operating elements

1	Adapter housing (PP)	7	Battery compartment
2	Adapter/non-return valve (SI/PTFE)	8	Operating mode selection switch
3	Membrane filter (PP/PTFE)	9	Power adapter: Input: AC 100 - 240 V; 50/60 Hz Output: DC 5V/200 mA
4	Pressure compensation opening	10	Plug
5	Filling button	11	Charging socket
6	Dispensing button		

4 First steps - Connecting the battery


- Grip the battery plug by the wire, confirm alignment of the plug and insert it into the socket. Push it firmly into the socket.

- b. Insert the battery.

NOTICE! A flashing LED indicates that the battery must be charged prior to use, see Charging the battery, p. 26.

- d. Close the battery compartment.

5 Pipetting

- a. Always hold the pipette as close as possible to the upper end and push it carefully into the adapter until it is firmly in place.

⚠ WARNING! Risk of injury!

Never use force! Especially with thin pipettes, there is a risk of glass breakage! Ensure that the pipette is firmly in place.

- b. Always hold the instrument with the pipette inserted vertically, tip down.

Filling


- a. Immerse pipette tip into the test liquid.
- b. Press the upper pipetting button slowly and fill the pipette so that the meniscus lies slightly above the desired marks.

NOTICE! Ensure that the pipette is not overfilled.

Changing the pipetting speed

- a. To fill, press the upper pipetting button; to dispense, press the lower pipetting button.

The pipetting speed increases continuously up to the set maximum motor speed the further the pipetting button is pressed in.

Setting the volume

- a. Wipe the tip of the pipette with a suitable lint-free cloth.

Press the lower pipetting button slowly and dispense liquid until the meniscus is adjusted exactly to the desired volume.

Selecting the operating mode

- a. Switch up: Free delivery
- b. Switch down: Blow-out (with motor power)

Emptying

- a. Press the lower pipetting button slowly. The further the pipetting button is pressed in, the faster the pipette is emptied.

Check tightness

- a. After adjusting the meniscus, no liquid should drip from the pipette without operating the valves. If the pipette drips; see Troubleshooting, p. 30

6 Charging the battery

NOTICE! Use only the original power adapter supplied by the manufacturer!

Using another power adapter can damage the instrument and the battery. The voltage specified on the power adapter must match the mains voltage. The mains voltage may deviate by a maximum of $\pm 10\%$ from the specification on the power adapter. Always verify this!

The built-in nickel-metal hydride battery guarantees an operating time of approx. 8 hours when fully charged. A flashing LED indicates that the battery needs to be recharged.

Charge indicator

LED off or flashing green:

sufficient battery capacity

LED flashes red slowly:

remaining battery capacity is sufficient for a max. of 2 hours of pipetting; charging the battery is recommended

LED flashes red quickly:

Stop the pipetting process and recharge the battery

Charging process

- Charging temperature: + 10 °C to + 35 °C
- After connecting the power adapter, the LED display of the instrument lights up continuously.
- After charging for four hours, the power adapter automatically switches to a pulsed trickle charging (LED display flashes every 15 seconds).

NOTICE! Trickle charging:

Trickle charging prevents overcharging and minimizes the lazy battery effect. To maintain maximum capacity, only recharge the battery when the LED charge indicator lights up. If the LED display does not light up during charging, see Troubleshooting, p. 30


6.1 Replacing the battery

⚠ WARNING! Incorrect batteries


Using incorrect batteries or handling the instrument improperly (short circuit, mechan. destruction, overheating, etc.) can lead to the battery exploding in extreme cases.


NOTICE! Use only original accessories, see Accessories and replacement parts, p. 31. Do not replace the battery with non-rechargeable batteries or rechargeable batteries of other manufacturers.

- a. To open the battery compartment, press on the arrow marking, position ① while sliding the cover backwards.


- b. grip the battery plug by the wire and gently pull it out of the socket.

- c. Grip the plug of the new battery by the wire, confirm alignment of the plug and push it firmly into the socket until you hear it snap into place.

- d. Insert a new battery and close the battery compartment.

7 Replace filter and clean adapter

⚠ WARNING

Risk of contamination and infection

Wear protective gloves and eye protection.


NOTICE! Filter replacement or cleaning is only necessary if a pipette has been overfilled or the suction capacity is impaired.

- a. Pull out pipette.
- b. Unscrew the adapter housing.
- c. Remove the wetted membrane filter.

- d. Press the adapter upwards out of the adapter housing. Do not remove the non-return valve from the adapter.

- e. For cleaning, rinse the adapter carefully with a wash bottle filled with a suitable cleaning liquid (e.g., water). Then blow out and dry.

- f. First press the new membrane filter into the adapter with the thicker connection facing downwards (the thinner, conical side facing upwards) and insert it into the adapter housing.
- g. The complete pipette holder incl. membrane filter are autoclavable at 121°C (2 bar) and a holding time of at least 15 minutes, in accordance with DIN EN 285.
- h. Reassemble the instrument in reverse order and check for leaks, see section Pipetting, p. 24.

NOTICE! Autoclavable filters:

The effectiveness of autoclaving must be verified by the user. Not every membrane filter is suitable! Use only original accessories.

- The instrument comes equipped with a 0.2 µm membrane filter (red color code).
- The membrane filter can be autoclaved up to five times at 121°C (250 °F).
- Class B volume tolerances can be achieved with 0.2 µm membrane filters (for cell cultures).

- 3 µm membrane filters are required to comply with Class A volume tolerances. Due to the larger pore size, these have a reduced impact on the free delivery time of the pipettes, see Accessories and replacement parts, p. 31

8 Cleaning | UV disinfection | Autoclaving

Occasionally clean the housing parts with a damp cloth.

When used properly, the instrument is maintenance-free.

The device is resistant to normal exposure to a UV disinfection lamp. The device may change color due to the UV exposure.

The membrane filter, adapter and adapter housing are autoclavable at 121 °C, 2 bar and a holding time of at least 15 minutes, in accordance with DIN EN 285.

Before autoclaving, the parts must be cleaned carefully, see Replace filter and clean adapter, p. 28.

9 Troubleshooting

Problem	Cause	Corrective action
Suction capacity impaired	• Non-return valve closed • Filter or non-return valve dirty	• Set the selection to “Blow-out” and press the lower pipetting button, see Pipetting, p. 24 • Clean the non-return valve and replace the filter if necessary (Replace filter and clean adapter, p. 28)
Pipette dripping	• Filter not properly positioned • Adapter or pipette damaged	• Insert the filter correctly • Replace the adapter or pipette, see Pipetting, p. 24 and Replace filter and clean adapter, p. 28, p. 28
Pipette not held in adapter	• Adapter dirty or damaged	• Clean adapter, allow to air dry or replace if necessary, see Replace filter and clean adapter, p. 28
After connecting the power adapter, the LED display of the instrument flashes or does not light up	• Battery/power adapter or instrument defective	• Replace battery/power adapter • Send instrument in for repair, see Sending for repair, p. 32
Charging process unsuccessful; LED dis-	• Battery defective	• Replace the battery, see Replacing the battery, p. 27

Problem	Cause	Corrective action
play of the instrument does not light up		

If the troubleshooting table above does not solve the problem, contact your authorized dealer or the manufacturer directly.

10 Ordering Information

pipeo®- pipette controllers with power adapter

Items supplied: VITLAB pipeo®, power adapter (100 - 240 V, 50/60 Hz), 4 interchangeable plugs (EU, UK, US/J, AUS), battery, one battery compartment cover, two replacement membrane filters (0.2 µm), and operating manual

Item no. 1631500

10.1 Accessories and replacement parts

Power adapter for Europe (continental) Item no. 1671831	Silicone adapter with non-return valve Item no. 1670203
Power adapter for UK/Ireland Item no. 1671832	Replacement membrane filter 0.2 µm, sterile (red color code) 1 pc in blister packaging Item no. 1670647
Power adapter for Japan Item no. 1671833	Replacement membrane filter 0.2 µm, non-sterile 10 pc in PE bag Item no. 1670648
Power adapter for Australia Item no. 1671834	Replacement membrane filter 3 µm, non-sterile 10 pc in PE bag Item no. 1670650
Power adapter for USA Item no. 1671835	Nickel-metal hydride battery Item no. 1670218
Adapter housing Item no. 1661281	Wall mount, gray Item no. 1670660

11 Warranty

We shall not be liable for the consequences of improper handling, use, servicing, operating or unauthorized repairs of the device or for the consequences of normal wear and tear, especially of wearing parts such as pistons, seals, valves and the breakage of glass. The same applies for failure to follow the instructions of the operating manual. We are not liable for damage resulting from disassembly beyond that described in the operating manual or if non-original spare parts or components have been installed.

12 Repair – calibration service

If any instrument malfunction cannot be remedied in the lab by simply exchanging replacement parts, the instrument must be sent in for repair.

Please note that for safety reasons only clean and decontaminated instruments can be inspected and repaired!

12.1 Sending for repair

NOTICE! Transporting hazardous materials without approval is prohibited by law.

1. Clean the instrument thoroughly and decontaminate.
2. Complete the “Declarations of Health Safety Compliance” form (forms can be requested from the dealer or manufacturer or downloaded from www.vitlab.com).
3. Send the completed form together with the instrument to the manufacturer or dealer. Include an exact description of the malfunction and the media used.

Shipment is at the risk and the cost of the sender.

13 Disposal

The symbol  means that at the end of their service life, batteries/accumulators and electronic devices must be disposed of separately from household waste (unsorted municipal waste).

Electronic devices must be disposed of in accordance with Directive 2012/19/ EU of the European Parliament and of the Council from July 04, 2012 on waste from electrical and electronic equipment and in compliance with national disposal regulations.

Both batteries and accumulators (rechargeable batteries) contain materials that can be damaging to the environment and human health. Therefore, they must be properly disposed of in accordance with Regulation (EU) 2023/1542 of the European Parliament and of the Council from July 12, 2023 on batteries and waste batteries and in compliance with national disposal regulations. Only dispose of fully discharged batteries and accumulators.

⚠ WARNING! Fire hazard due to heating after a short circuit

Do not short-circuit the battery to discharge it! Tape the plug with adhesive tape to prevent a short circuit in the collection container! Never disassemble the battery!

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