Pipeta / Pipette

VITLAB®

micropipette

Instrucciones de manejo
Mode d'emploi

VITLAB GmbH
Linus-Pauling-Str.1
63762 Grossostheim
Germany
tel: +49 6026 97799-0
fax: +49 6026 97799-30
info@vitlab.com
www.vitlab.com

Contenido

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
7. Pipetting


- 5 ml and 10 ml instruments should only be used with the PE filter installed (see page 30).
- Pipette tips are disposable items!

1. Fitting the tip


Use the correct tips according to the volume range or the color code.

Ensure that the tip is securely seated.

2. Volume setting

Select the desired volume by rotating the volumes setting wheel. Avoid twisting and abrupt rotating motions during this adjustment.

3. Aspirate sample


a) Press pipetting button to the first stop.


b) Hold the pipette vertically and immerse the tip 2 to 3 mm into the liquid.

| Volume range | Immersion depth in mm | Waiting time in s |
|------------------------------|-----------------------|-------------------|
| > 1 μ l - 100 μ l | 2 - 3 | 1 |
| > 100 μ l - 1000 μ l | 2 - 4 | 1 |
| > 1000 μ l | 3 - 6 | 3 |

c) Let the pipetting button slide back slowly. In order for the liquid to reach its end position, leave the tip immersed for another 1-2 s.

d) Touch the tip against the container wall.

4. Discharge sample


- Place the pipette tip against the wall of the vessel. Hold the pipette at an angle of 30-45° relative to the container wall.
- Press the pipetting button slowly to the first stop and hold it down. For serum and liquids of high viscosity or low surface tension, observe adequate waiting time to improve accuracy.

- The blow-out stroke empties the tip completely: Press the pipetting button down to the second stop.
- While doing this, wipe the pipette tip against the wall of the container.
- Remove the pipette tip from the container wall and let the pipetting button slide back.

5. Ejecting the tip

Hold the pipette shaft over a suitable disposal container and press the tip ejection key to the stop.

Note:

ISO 8655 prescribes rinsing the pipette tip once with the sample liquid prior to the actual pipetting process.

Important!

Don't lay the instrument horizontal when the tip is filled. Liquid may enter and contaminate the instrument. The instrument should be stored without tips, placed upright in the shelf/rack mount or bench top rack which can be ordered separately.

8. Checking the Volume

Depending on use, we recommend inspection of the instrument every 3 to 12 months. The cycle can, however, be adjusted to individual requirements.

The gravimetric testing of the pipette volume is performed according to the following steps and is in accordance with DIN EN ISO 8655, Part 6.

1. Set nominal volume

Set volume to the maximum volume indicated on the instrument (see page 25 for procedure).

2. Condition the pipette

Condition the pipette before testing by using a pipette tip to aspirate and discharge the test liquid (distilled H₂O) five times.

3. Carry out the test

- Attach new pipette tip and prerinse one time with test liquid.
- Aspirate liquid and pipette it into the weighing vessel.

- Weigh the pipetted quantity with an analytical balance. Please follow the operating manual instructions from the balance manufacturer.
- Calculate the volume, taking the temperature of test liquid into account.
- At least 10 pipettings and weighings in three volume ranges (100%, 50%, 10% of nominal volume) are recommended for statistical analysis.

Calculation for nominal volume V_0

x_i = Weighing results

n = Number of weighings

Z = Correction factor (e.g. 1,0029 µl/mg at 20 °C, 1013 hPa)

Mean value

$$\bar{x} = \frac{\sum x_i}{n}$$

Mean volume

$$\bar{V} = \bar{x} \cdot Z$$

Standard Deviation

$$s = Z \cdot \sqrt{\frac{\sum (x_i - \bar{x})^2}{n - 1}}$$

Accuracy*

$$A\% = \frac{\bar{V} - V_0}{V_0} \cdot 100$$

Coefficient of Variation*

$$CV\% = \frac{100 \cdot s}{\bar{V}}$$

*) = Calculation of accuracy (A%) and variation coefficient (CV%): A% and CV% are calculated according to the formulas for statistical control.

Note:

Testing instructions (SOPs) are available for download at www.vitlab.com.

9. Accuracy Table

VITLAB® micropipette, Digital adjustable

| Volume range µl | Volume step µl | A* ≤ ± % | CV* ≤ % | Increment µl | Recommended type of tip, µl |
|--------------------|-------------------|-------------|------------|-----------------|--------------------------------|
| 0.5 - 10 | 10 | 1 | 0.5 | 0.01 | 0.5 - 20 |
| | 5 | 1.6 | 1 | | |
| | 1 | 7 | 4 | | |
| 2 - 20 | 20 | 0.8 | 0.4 | 0.02 | 2 - 200 |
| | 10 | 1.2 | 0.7 | | |
| | 2 | 5 | 2 | | |
| 10 - 100 | 100 | 0.6 | 0.2 | 0.1 | 2 - 200 |
| | 50 | 0.8 | 0.4 | | |
| | 10 | 3 | 1 | | |
| 20 - 200 | 200 | 0.6 | 0.2 | 0.2 | 2 - 200 |
| | 100 | 0.8 | 0.3 | | |
| | 20 | 3 | 0.6 | | |
| 100 - 1000 | 1000 | 0.6 | 0.2 | 1 | 50 - 1000 |
| | 500 | 0.8 | 0.3 | | |
| | 100 | 3 | 0.6 | | |
| 500 - 5000 | 5000 | 0.6 | 0.2 | 5 | 500 - 5000 |
| | 2500 | 0.8 | 0.3 | | |
| | 500 | 3 | 0.6 | | |
| 1000 - 10000 | 10000 | 0.6 | 0.2 | 10 | 1000 - 10000 |
| | 5000 | 0.8 | 0.3 | | |
| | 1000 | 3 | 0.6 | | |

Final test values related to the nominal capacity (maximum volume) indicated on the instrument, obtained when instrument and distilled water are equilibrated at ambient temperature (20 °C/68 °F) and with smooth operation. According to DIN EN ISO 8655.

* A = Accuracy

* CV = Coefficient of Variation

DE-M


20 °C
Ex

10. Adjustment


The instrument is permanently adjusted for aqueous solutions. If the pipette operation is clearly inaccurate, or if the instrument must be adjusted for solutions of different densities and viscosities or specially-shaped pipette tips, adjustments can be made.

1. Check the volume, determine actual value (see page 27).


2. Remove the cover: Push the hook forward, raise it slightly and then pull it back.


3. Using a paperclip or an unused pipette tip, remove the protective film (this protective film can be discarded).

4. Push the red adjustment slider completely back, raise the volume-setting wheel (decoupling) and release the adjustment slider.

5. Set the adjustment value:
set the volume-setting wheel to the previously determined actual value. A volume check is recommended after every adjustment.

6. Push the adjustment slider completely back again, push the volume-setting wheel downwards and release the adjustment slider. Re-insert the cover.

Note:

The change to the factory settings is indicated by the red adjustment slider now visible in the label window.

11. Autoclaving

The pipette is completely autoclavable at 121 °C (250 °F), 2 bar absolute (30 psi) with a holding time of at least 15 minutes, according to DIN EN 285.

1. Eject the pipette tip.
2. Autoclave the complete pipette without any further disassembling. Remove the filter of pipette 5 ml and 10 ml before autoclaving.
3. Allow the pipette to completely cool and dry.

Note:

The effectiveness of the autoclaving must be verified by the user. Maximum reliability is obtained with vacuum sterilization. We recommend the use of sterilization bags.

Attention!

Prior to autoclaving, the volume adjustment must be set on an available numbered volume (e.g., 11.25 or 11.26 but not between).

If the pipette is autoclaved frequently, the piston and the seal should be greased with silicone grease in order to preserve smooth movement. Please use only the recommended silicone grease, see accessories page 33. If necessary after sterilization, tighten the connection between the hand grip and the pipette shaft.

12. UV sterilization

The unit can withstand the usual output of a UV sterilization lamp. The effects of the UV may cause some color change.

13. Filter pipette 5 ml / 10 ml

A hydrophobic PE filter is used as a safeguard against liquid entering the pipette.

Change the filter if it becomes wet or contaminated.

- Use a flat object such as a screwdriver
- Remove the filter without damaging the tip cone.

Remove the filter before autoclaving!

The instrument can be operated without a filter.

14. Servicing and Cleaning

14.1. VITLAB® micropipette up to 1000 µl

1. Servicing

Inspect the pipette tip cone for damage.

Inspect the piston and seal for contamination.

Test the instrument's piston seal. To do this affix a tip, and aspirate a sample. Hold the instrument vertically, with the sample in the tip for approximately 10 s. If a drop forms at the tip orifice, see the troubleshooting guide on page 36.


2. Disassembly and cleaning

1. Unscrew the pipette shaft (S) from the hand grip.
2. Unscrew the upper part of the ejector (A) from the pipette shaft.
3. Pull the shaft (B, C and D) out of the lower part (E) of the ejector.
4. Unscrew the piston unit (B).

Note:

Piston remains connected with piston unit (B)!

5. Remove the seal with spring (C) (this is non-removable on 10 µl pipette models).
6. Clean the parts shown with a mild soap solution or isopropanol and then rinse with distilled water.
7. Allow the parts to dry (max. 120 °C/ 248 °F).
8. Grease piston and seal with a very thin layer of silicone grease. Please use only the recommended silicone grease, see accessories page 33.
9. Assemble the ambient temperature parts in reverse order from above. Piston unit and upper part of the ejector (A, B) should only be hand-tight.

(For illustration purposes only)

14.2. VITLAB® micropipette 5 ml and 10 ml

1. Servicing

Inspect the pipette tip cone for damage.

Inspect the piston and O-Ring-seal for contamination.

Test the instrument's piston seal. To do this, affix a tip, and aspirate a sample. Hold the instrument vertically, with the sample in the tip for approximately 10 s. If a drop forms at the tip orifice, see the troubleshooting guide on page 36.

2. Disassembly and cleaning

1. Remove the entire shaft (S) from the handgrip by rotating at the upper end of the ejector (F) and remove the filter (K) from the bottom part of the shaft (H).
2. Separate the bottom part of the ejector (F') by unscrewing it from the upper part of the ejector (F).
3. Unscrew and dismantle the piston unit (G) with the ejector spring (I) and the bottom part of the shaft (H).
4. Remove the O-Ring-seal from the piston unit and clean it.

Note:

Do not disassemble piston unit (G) any further!

5. Clean piston unit (G) and lower part of pipette shaft (H) with a soap solution or isopropanol and then rinse with distilled water.
6. Allow the parts to dry (max. 120 °C / 248 °F) and to cool down.
7. Carefully lubricate the inside and outside of the O-ring and mount it on the piston.
8. Assemble the individual components in the reverse order from that described above.

(For illustration purposes only)

15. Ordering Information and Accessories

15.1. Ordering Information

VITLAB® micropipette

| Volume | Cat. No. |
|---------------|----------|
| 0.5 - 10 µl | 1641000 |
| 2 - 20 µl | 1641002 |
| 10 - 100 µl | 1641004 |
| 20 - 200 µl | 1641006 |
| 100 - 1000 µl | 1641008 |
| 0.5 - 5 ml | 1641010 |
| 1.0 - 10 ml | 1641012 |


15.2. Accessories

Bench-top rack for 6 pipettes

Cat. No. 1672002

Shelf/rack mount for 1 pipette

Cat. No. 1672000

Filter for VITLAB® micropipette 5 ml,

pack of 25

Cat. No.1672010

Filter for VITLAB® micropipette 10 ml,

pack of 25

Cat. No. 1672012

Silicone grease for VITLAB® micropipette up to 1000 µl

with stainless steel piston, pack of 1 Cat. No. 1672015

Silicone grease for VITLAB® micropipette 5 ml/10 ml

pack of 1


Cat. No. 1672016

16. Spare parts

16.1. VITLAB® micropipette up to 1000 µl

Before reordering a piston unit or seal and spring for your VITLAB® micropipette with glass piston (up to serial number 08N), please contact info@vitlab.com.

(For illustration purposes only, parts will differ slightly depending on nominal volume of instrument.)


| Volume | A | B | C | D | E |
|---------------|---------|---------|---------|----------|---------|
| 0.5 - 10 µl | 1671400 | 1671411 | - | 1671431* | 1671441 |
| 2 - 20 µl | 1671401 | 1671412 | 1671420 | 1671432 | 1671442 |
| 10 - 100 µl | 1671402 | 1671416 | 1671424 | 1671433 | 1671443 |
| 20 - 200 µl | 1671403 | 1671417 | 1671425 | 1671434 | 1671444 |
| 100 - 1000 µl | 1671404 | 1671418 | 1671426 | 1671435 | 1671445 |

* 0.5-10 µl including seal

16.2. VITLAB® micropipette 5 ml and 10 ml

(For illustration purposes only, parts will differ slightly depending on nominal volume of instrument.)

Upper part
of ejector

Lower part
of ejector

Piston
unit

Lower part
of pipette
shaft

Ejector
spring

| Volume | F + F' | G | H | I |
|------------|---------|---------|---------|---------|
| 0.5 - 5 ml | 1671451 | 1671455 | 1671461 | 1671465 |
| 1 - 10 ml | 1671452 | 1671456 | 1671462 | 1671465 |

17. Troubleshooting

| Problem | Possible cause | Corrective action |
|--|---|--|
| Tip dripping (instrument leaks) | <ul style="list-style-type: none">- Unsuitable tip- Tip not seated tightly | <ul style="list-style-type: none">- Only use high-quality tips- Press tip on firmly |
| The instrument does not aspirate or aspirates too little; the discharged volume is too low | <ul style="list-style-type: none">- Seal contaminated- The seal or cone is damaged- The piston is contaminated or damaged | <ul style="list-style-type: none">- Clean seal- Replace seal or shaft- Clean or replace piston |
| Aspiration is too slow | <ul style="list-style-type: none">- Shaft clogged- The filter in the 5 ml and 10 ml models is contaminated | <ul style="list-style-type: none">- Clean shaft- Change the filter |
| Discharged volume is too large | <ul style="list-style-type: none">- Pipetting button pressed too far into the blow-out position before sample uptake | <ul style="list-style-type: none">- Operate properly. See 'Pipetting', page 25. |
| Piston is difficult to move | <ul style="list-style-type: none">- The piston is contaminated or needs oil | <ul style="list-style-type: none">- Clean and oil piston |

18. Repairs - Calibration Service

If a problem cannot be fixed by following the troubleshooting guide, or by replacing spare parts, then the instrument must be sent in for repair.

For safety reasons, instruments returned for checks and repairs must be clean and decontaminated!

18.1. Return for repair

- a) Clean and decontaminate the instrument carefully.
- b) Complete the „Declaration on Absence of Health Hazards“ (ask your supplier or manufacturer for the form. The form can also be downloaded from www.vitlab.com).
- c) Send the completed form along with the instrument to the manufacturer or to the dealer with an exact description of the type of malfunction and the media used.

The return transport of the instrument is at risk and cost of the sender.

18.2. Calibration Service

ISO 9001 and GLP-guidelines require regular examinations of your volumetric instruments. We recommend checking the volume every 3-12 months. The interval depends on the specific requirements on the instrument. For instruments frequently used or in use with aggressive media, the interval should be shorter. The detailed testing instruction can be downloaded on www.vitlab.com.

VITLAB also offers you the possibility to have your instruments calibrated by the VITLAB Calibration Service.

Just send in the instruments to be calibrated, accompanied by an indication of which kind of calibration you wish. Your instruments will be returned within a few days together with a test report (VITLAB calibration service) or with a DAkkS Calibration Certificate. For further information, please contact your dealer or VITLAB.

19. Warranty

We shall not be liable for the consequences of improper handling, use, servicing, operating or unauthorized repairs of the instrument or the consequences of normal wear and tear especially of wearing parts such as pistons, seals, valves and the breakage of glass as well as the failure to follow the instructions of the operating manual. We are not liable for damage resulting from any actions not described in the operating manual or if non-original parts have been used.

20. Disposal

For the disposal of instruments and tips, please observe the relevant national disposal regulations.

Subject to technical modification without notice.
We will not be held responsible for printing or typographical errors.